

**Center for Veterinary Biologics
and
National Veterinary Services Laboratories
Testing Protocol**

**Supplemental Assay Method for the Determination of
Protein and Phenol in PPD (Purified Protein Derivative
Produced From Cultures of *Mycobacterium bovis* Strain
AN-5) Bovis Tuberculin**

Date: July 26, 2000
Supersedes: December 10, 1998
Number: TCSAM0513.05
Standard Requirement: 9 CFR, Part 113.409
Contact Person: Gary B. Clark, (515) 663-7818

Approvals:

P. Frank Ross, Head
Toxicology and Chemistry Section
Date: _____

Arthur J. Davis, Chief
Pathobiology Laboratory
Date: _____

Ann L. Wieggers, Quality Assurance Manager
Date: _____

Randall L. Levings, Director
Center for Veterinary Biologics-Laboratory
Date: _____

/s/ William Buisch _____ Date: 8/16/00__
William A. Buisch, Acting Director
National Veterinary Services Laboratories

United States Department of Agriculture
Animal and Plant Health Inspection Service
P. O. Box 844
Ames, IA 50010

Mention of trademark or proprietary product does not constitute a guarantee or warranty of the product by USDA and does not imply its approval to the exclusion of other products that may be suitable.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

Table of Contents

1. Introduction
 - 1.1 Background
 - 1.2 Key words
2. Materials
 - 2.1 Equipment
 - 2.2 Reagents/supplies
 - 2.2.1 Protein test
 - 2.2.2 Phenol test
3. Preparation for the test
 - 3.1 Personnel qualifications/training
 - 3.2 Preparation of equipment/instrumentation
 - 3.3 Preparation of reagents/control procedures
 - 3.3.1 Protein test
 - 3.3.2 Phenol test
 - 3.4 Preparation of the sample
 - 3.4.1 Receipt
 - 3.4.2 Preparation
4. Performance of the test
 - 4.1 Protein
 - 4.2 Phenol
5. Interpretation of the test results
 - 5.1 Protein
 - 5.2 Phenol
 - 5.3 Controls
6. Report of test results
7. References
8. Summary of revisions

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

1. Introduction

1.1 Background

The Code of Federal Regulations, Title 9 (9 CFR) (Animals and Animal Products) states that the Animal and Plant Health Inspection Service (APHIS) is responsible for administering the Virus-Serum-Toxin Act. It specifies testing methods for licensed tuberculin products. Protein concentration is determined by classical Kjeldhal digestion, distillation, and titration of the ammonia. Phenol is determined by end-point titration with bromate/bromide. Satisfactory product must contain $1.0 \text{ mg/ml} \pm 0.1 \text{ mg/ml}$ protein. Phenol content must be $0.50\% \pm 0.04\%$.

1.2 Key words

protein, phenol, PPD, tuberculin

2. Materials

2.1 Equipment

2.1.1 Balance, top loading, capable of measuring 0.01 g

2.1.2 Digestion unit, Buchi, B-426, with digestion tubes, or equivalent

2.1.3 Distillation unit, Buchi, B-316, or equivalent

2.1.4 Volumetric pipets, Class A, meets ASTM Standard E969-83

2.1.5 Volumetric flasks, Class A, with barrel head glass stopper, meets ASTM E288 requirements

2.1.6 Erlenmeyer flasks, 125 ml

2.1.7 Buret with PTFE stopcock, 10 ml, precision bore, calibrated to ASTM E-694 accuracy requirements

**Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin**

2.1.8 Buret with PTFE stopcock, 50 ml, precision bore, calibrated to ASTM E-694 requirements

2.1.9 Graduated cylinders, 50, 100, 250, 500, and 1,000 ml (PYREX), meeting ASTM D86, D216, and D447 requirements

2.1.10 Glass-stoppered erlenmeyer flasks, 250 ml

2.1.11 Heating/stirring plate with stirring bars

2.1.12 Fast filter paper, Whatman No. 1

2.2 Reagents/supplies

All chemicals, reagent grade. Use distilled or demineralized water or water of equivalent purity.

2.2.1 Protein test

- 1.** Sulfuric acid (H_2SO_4)--Purity: Minimum 95.0%, Maximum 98.0%
- 2.** Mercury tablets, Brinkmann Instruments, Catalog No. 015-00-646-3
- 3.** Sodium hydroxide (NaOH)--Purity: 98.5%
- 4.** Boric acid (H_3BO_3)--Purity: 99.9%
- 5.** Methyl red--Purity: 98.0%
- 6.** Hydrochloric acid (HCl)--Assay: 36.5%-38.0%
- 7.** Sodium carbonate (Na_2CO_3)--Purity: 99.9%
- 8.** Bromo phenol blue--Purity: 98.0%
- 9.** National Veterinary Services Laboratories (NVSL) Control--Pool of PPD tuberculin products with established protein and phenol values

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

10. Protein-Standard, National Institute of Standards and Technology, Gaithersburg, MD 20899, Standard Reference Material® 927 C, Bovine Serum Albumin, Certified Protein Concentration 71.57 g/L

2.2.2 Phenol test (some reagents same as for protein)

1. Methyl orange--Purity: 98.0%
2. Silicotungstic acid ($\text{H}_4[\text{Si}(\text{W}_3\text{O}_{10})_4] \cdot 26\text{H}_2\text{O}$)--Purity: 99.0% Store at 4°C.
3. Arsenic trioxide (As_2O_3)--Purity: 99.9%
4. Sodium bicarbonate (NaHCO_3)--Purity: 99.9%
5. Potassium bromate (KBrO_3)--Purity: 98.5%
6. Potassium bromide (KBr)--Purity: 99.0%
7. Phenol ($\text{C}_6\text{H}_5\text{OH}$)--Purity: $\geq 99.0\%$

3. Preparation for the test

3.1 Personnel qualifications/training

No special test-related training is needed for this testing. Analysts performing this procedure should first conduct 2 trial runs using controls and standards and obtain results within acceptable limits.

3.2 Preparation of equipment/instrumentation

Become familiar with Buchi instruction regarding operation. Turn on water that aspirates fumes from the digestion unit and keeps the water cool in the condenser of the distillation unit. Adjust water flow in the distillation unit to approximately 1 L per min. Turn on the distillation unit. Set time preselector to "2" (2 min) and stopcock for aspiration to "Off." Make sure that Buchi bottles of NaOH and water are adequately filled.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

3.3 Preparation of reagents/control procedures

3.3.1 Protein test (all reagents stable for at least 6 mo unless specified)

1. Cut Hg tablets in half.

Caution: Because tablets contain mercury, handle in fume hood and wear gloves, protective glasses, and mask.

2. 32% NaOH, dissolve 640 g \pm 1 g NaOH in 1.4 L H₂O in 2-L volumetric flask on the magnetic stirrer. Cool to room temperature (RT). Dilute to volume with H₂O. Repeat above until Buchi 10-L bottle is full. Store at RT.

Caution: NaOH is caustic--Avoid contact with skin.

3. Saturated H₃BO₃, add 15 g to 100 ml H₂O. Stir, with heat, until all H₃BO₃ dissolves. Some H₃BO₃ recrystallizes when cool. Store at RT.

4. 0.1% bromo phenol blue, dissolve 0.1 g in 100 ml H₂O. Store at RT.

5. 0.5% methyl red, dissolve 0.5 g in 100 ml ethanol. Store at 4°C.

6. Standardized 0.01 N HCl-0.02 N HCl, 1.7 ml HCl/L H₂O. Titrate approximately 0.0100 g dried sodium carbonate dissolved in 25 ml H₂O. Indicator: 3 drops 0.1% bromo phenol blue; the color of endpoint is green, not bluish green nor yellowish green. Store at RT.

Calculation:

$$N \text{ HCl} = [(g \text{ Na}_2\text{CO}_3) \times (1000)] / [(Vol \text{ HCl}) \times (52.994)].$$

Caution: Concentrated HCl is corrosive--Handle in fume hood. Avoid contact with skin.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

7. Protein Standard, dilute protein (Section 2.2.1.10) to the range of 0.9-1.1 mg/ml. Prepare sufficient dilution to provide several aliquots of 15-ml portions in 30-ml serum vials. Store at 4°C.

3.3.2 Phenol test (all reagents stable for at least 6 mo unless specified)

1. 20% HCl, slowly add 200 ml HCl to 600 ml H₂O; dilute to 1 L. Store at RT.
2. 0.1% methyl orange, add 0.1 g methyl orange to 100 ml H₂O. Filter if necessary. Store at RT.
3. Silicotungstic acid solution (SAS), dissolve 60 g H₄[Si(W₃O₁₀)₄]*26H₂O in 400 ml H₂O in 500-ml volumetric flask. Add 50 ml H₂SO₄. When cool, dilute to volume with H₂O. Store at RT.
4. Clarifying solution (CS), add 50 ml SAS and 125 ml 20% HCl to 325 ml H₂O. Prepare fresh prior to each test.
5. "Acid solution" for As₂O₃ standard solution, add 110 ml HCl and 2.5 ml methyl orange to 100 ml H₂O. Store at RT.
6. 0.0500 N As₂O₃, dissolve 2.4730 g dried As₂O₃ in 25 ml hot 1N NaOH in 1-L volumetric flask. Neutralize with 25 ml 1N H₂SO₄. Cool and dilute to volume with H₂O. Store at RT.

Caution: As₂O₃ is extremely toxic--Avoid contact; handle in fume hood using gloves, mask, and goggles. Consult Material Safety Data Sheet for specific handling instructions.

7. Phenol standard, 0.50%, dissolve 2.50 g phenol in 500 ml H₂O. Prepare in 500-ml volumetric flask. Store at RT.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

Critical Control Point: The final diluted volume of the test fluid must be adjusted as described in Section 4.3.2.8.

8. Test fluid (TF), dissolve 0.30 g NaHCO_3 , 1.67 g KBrO_3 , and 15.00 g KBr in H_2O and qs to 1 L with H_2O . Store at RT. The TF volume must be adjusted by adding corrected volume of H_2O to TF. It must take a volume of 21.3 ml to titrate 25 ml 0.050 N As_2O_3 in 10 ml "Acid Solution." A first time titration will require less than 21.3 ml TF. Adjust as described in the following example:

Example: Assume the first time titration volume is 20.5 ml.

$$(1,000 \text{ ml of TF}) - (20.5 \text{ ml}) = 979.5 \text{ ml}$$

$$\frac{(979.5)(\text{desired vol})}{(\text{actual vol})} \text{ or } \frac{(979.5)(21.3)}{(20.5)} = 1,017.2 \text{ ml}$$

For corrected volume of H_2O : $1017.2 - 979.4 = 37.8 \text{ ml}$ to be added to TF.

Note: TF in buret has to be put back into flask.

3.4 Preparation of the sample

3.4.1 Receipt

Reference current version of **TCSOP0001**.

3.4.2 Preparation

PPD tuberculin products are received in sealed serum bottles. They are stored at 4°C in the walk-in refrigerator prior to testing. Before testing, allow sample vials and reagents to warm to room temperature.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

4. Performance of the test (use current version of TCFRM0513)

4.1 Protein

(Analyze the control pool and protein standard each time testing is performed. Analyze each in triplicate.)

4.1.1 Place 5.0 ml sample, one half Hg tablet, and 3.0 ml H₂SO₄ into a digestion tube

Caution: HgO is poisonous--Use gloves, mask, and goggles.

Caution: Concentrated H₂SO₄ is corrosive--Avoid contact with skin.

4.1.2 Place the tubes in a digestion tube holder. Place the holder into the digestion unit. Turn on the unit and set energy regulator to "5." Fifteen min later, set to "7."

4.1.3 Digest until acid comes to true boil or no longer "burned smoke," about 50-60 min. Set to "9" for 15 more min.

4.1.4 Cool, add 6 ml H₂O, mix, and cool again.

4.1.5 Place digestion tube and a flask containing 5 ml H₃BO₃ and 3 drops indicator into the distillation unit. Tilt the flask so the tip of the condenser is immersed in the H₃BO₃.

4.1.6 Press and hold NaOH button and count to 3. Then hit Start button to start distillation unit. Distill for 2 min.

4.1.7 Titrate collected distillate to endpoint color change of yellow to deep rose (pH 5.0) with standardized HCl. Record the volume of HCl on log sheet.

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

4.2 Phenol

(Analyze the control pool and phenol standard each time testing is performed. Analyze each in triplicate.)

4.2.1 Add 5 ml sample and 100 ml CS to 250-ml glass-stoppered flask. Shake 2 min. Filter through filter paper into 50-ml cylinder.

4.2.2 Transfer 50 ml filtrate to another flask. Add 1 drop methyl orange, stopper, and shake a few sec. Observe the color; when red, go to **Section 4.2.3**.

4.2.3 Titrate with 2 ml test fluid (TF), stopper, and shake a few sec and observe the color. When red, repeat **Section 4.2.3**. When colorless, go to **Section 4.2.4**.

4.2.4 Shake 30 sec. Add 1 drop indicator, stopper, and shake a few sec and observe the color. When it does not turn to colorless within 10 sec, titrate with 1 ml TF, stopper, and repeat **Section 4.2.4**. When colorless, go to **Section 4.2.5**.

4.2.5 Shake 1 min. Add 1 drop indicator, stopper, and shake a few sec. Observe the color. When red stays longer than 10 sec, titrate with 0.50 ml TF, stopper, and repeat **Section 4.2.5**. When colorless, record total volume of TF as the endpoint of titration and use for calculation of percent phenol.

5. Interpretation of the test results

5.1 Protein (Report average of triplicates.)

$\text{mg Protein/ml} = (\text{ml HCl})(\text{N HCl})(1.400)(6.25)/(5 \text{ ml PPD})$

Satisfactory Protein Content: 1.0 mg/ml \pm 0.1 mg/ml

5.2 Phenol (Report average of triplicates.)

$\text{Percent phenol} = (\text{vol of test fluid})(0.04)-(0.04)$

Satisfactory Phenol Content: 0.50% \pm 0.04%

Supplemental Assay Method for the Determination of Protein and Phenol in PPD
(Purified Protein Derivative Produced From Cultures of *Mycobacterium bovis*
Strain AN-5) Bovis Tuberculin

5.3 Controls

Results for controls and standards must be within acceptable limits; otherwise repeat testing.

6. Report of test results

Validate and report results according to the current version of TCSOP0001.

7. References

7.1 Code of Federal Regulations, Title 9, Part 113.409, Revised as of January 01, 2000, page 653.

7.2 Official Methods of Analysis of AOAC International, Arlington, Virginia, 16th Edition, Pat Cuniff, Editor (1995), Volume I, Chapter 12, page 7.

8. Summary of revisions

8.1 Version .01 was written to meet NVSL/CVB Quality Assurance requirements, to clarify practices in use in the NVSL/CVB-L, and to provide additional detail. No significant changes were made from the previous protocol.

8.2 Version .02 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail.

8.3 Version .03 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail. The following are the significant changes made from the superseded protocol:

Change in the digestion apparatus

8.4 Version .04 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail.

8.5 Version .05 was written to clarify practices in use in the NVSL/CVB-L and to provide additional detail.